

PSMG2-proteasome (prosome, macropain) assembly chaperone 2 Gene

Reference: PTXBC013356

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CLAST3; HCCA3; HsT1707; MDS003; PAC2; TNFSF5IP1; proteasome assembly chaperone 2; CD40 ligand-activated specific transcript 3; HDCMC29P; HSPC260; PAC-2; hepatocellular carcinoma susceptibility protein; hepatocellular carcinoma-susceptibility protein 3; li		

For research use only.

GENE DESCRIPTION	proteasome (prosome, macropain) assembly chaperone 2
GENE ID	56984
ORIGIN SPECIES	Human
ACCESSIONS	BC013356
VECTOR	pDONR223
NCBI SYMBOL	PSMG2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcctccctgcgggagtcggcccccacattgccggtccacccctctaatgccagcagtcgtgtggaaatgttgccagcttgcaatggatctgattattctacactgaatatgtctaagattggttactctataccgattgtcttgccaatggttgaaacaatccatagcgacca
cagaaggaaatcaacagaacttagcataaatgctgaagtgatcattgccttcaagaaagctgggtgctctacagtttaagatccattttatagataaaatcaaacccgttctgtgaaaaactgcttccctgggtgaaaagcagtggtgctgccaagagtcattgttctcaagcagtcattcata
tcagcgtaatgatctgcagcttcgtagctacccctccggtacctaactacacctccatgcaaaaaagtgttcaaaataaaataaagagcccttaactgggaagaatgaaaaaacccgggtgattcctgaaatgatgatccgagtttgcacgcattccgggagaggtatcacaaaaac
actctatgatgaagctgttcaagaatccaaatggcagtcgtcgaattttgtcagaaggggacaacatccagatgcattaggtctgttgatctctaatgagtggttcagatactcaaacacttagcgcagacccacagatctgcctcacggtggaaaataccaagttcttga
gattactcttggcagtggttccccctgcattttctga