

CALU-calumenin Gene

Reference: PTXBC013383

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IEF SSP 9302; crocalbin-like protein; multiple EF-hand protein		
For research use only.			

GENE DESCRIPTION	calumenin
GENE ID	813
ORIGIN SPECIES	Human
ACCESSIONS	BC013383
VECTOR	pDONR223
NCBI SYMBOL	CALU
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacctgcagagtttcttatgtgcctgtccctgtgcacagcctttgcttgagcaaacccacagaaaaagagacgtgtacatcatgagcctcagctcagtgacaagggtcacaatgatgctcagagtttgattatgacctatgccttctgggtgctgaagaagcaaacgctttgat
cagctgacaccagagaagagcaaggaaagccttggaagattgtaagtaaaatagatggcgacaaggacgggtttgtcactgtggtgagctcaaaagactggalttaattgcacaaaagcgtggatttacgaggatgtagagcgacagtggaaagggcatgacctcaatgaggacg
gcctcgtttcctgggaggagtataaaatgccacctacggctacgttttagatgatccagatcctgatgatggatttaactataaacagatgatggttagatgagcggaggtttaaatggcagacaaggatggagacctatgccaccaaggaggatgttcacagcttctgcacctga
ggagtatgactacatgaaagatatagtagtacaggaacaatggaagatagataagaatgctgatggtttcattgatctagaagagtatatgttgatgacatgtacagccatgatgggaatactgatgagccagaatgggtaaagacagagcgagagcagttttgtgagtttcgggataagaa
ccgtgatgggaagatggacaaggaagaccaaagactggatccttccctcagactatgatcatgcagaggcagaagccagggcacctgtctatgaatcagacaaaacaaggatggcaagcttaccaaggaggagatcgttgacaagtatgactattttgtggcagccagccacag
attttgggagcgcttagtgcgcatgatgattctga