

CLEC7A-C-type lectin domain family 7, member A Gene

Reference: PTXBC013385

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BGR; CANDF4; CD369; CLECSF12; DECTIN1; SCARE2; C-type lectin domain family 7 member A; C-type (calcium dependent, carbohydrate-recognition domain) lectin, superfamily member 12; C-type lectin superfamily member 12; DC-associated C-type lectin 1; beta-gluc		

For research use only.

GENE DESCRIPTION	C-type lectin domain family 7, member A
GENE ID	64581
ORIGIN SPECIES	Human
ACCESSIONS	BC013385
VECTOR	pDONR223
NCBI SYMBOL	CLEC7A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaatcatcctgatttagaaaattggatgaagatggatatactcaattacacttcgactctcaagcaataaccaggatagctgtgttcagagaaaggatcgtgtgctgcactcctcttgccgcctcattgctgtaattttgggaatcctatgcttggaatactggatagctgtgtcct
gggtaccatggctgtgttcaaaagctgtgggaattcaaaaggataa