

RASD2-RASD family, member 2 Gene

Reference: PTXBC013419

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MGC:4834; TEM2; GTP-binding protein Rhes; Ras homolog enriched in striatum; tumor endothelial marker 2; RASD family member 2		

For research use only.

GENE DESCRIPTION	RASD family, member 2
GENE ID	23551
ORIGIN SPECIES	Human
ACCESSIONS	BC013419
VECTOR	pDONR223
NCBI SYMBOL	RASD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgaagactttgtccagcgggaactgcacgtcagtggtccgcaaaaactcaccgcatggtggtgctgggtgcctctcgggtggcgaagctccatcgtgtctcgtctcaatggcgccttgaggaccagtacacaccaccatcaggacttcaccgtaaggtatacaacatccgcggcgacatgtaccagctcgacatcctggatacctctggcaaccacccttccccgccatgcgcaggctgtccatcctcacaaggatgtcttcacctcgtgttcagcctggataaccggagtccttcgatagggtcaagccttcagaagcagatcctggaggtcaagtcctgcctgaagaacaagaccaaggaggcggcggagctgccccatgtcctctgtggcaacaagaacgaccacggcgagctgtgcccgccagggtgccaccaccaggccgagctgctgtgtcgggcgacgagaactgcgcctacttcgaggtgtcggccaagaagaacaccaacgtggacgagatgttctacgtccttcagcatggccaagctgccacacgagatgagcccgccctgcatcgcaagatctccgtgcagtacgtgacgccttcacccaggccctctcgtatgcgcgcgtcaaggagatggacgcctatggcatgtctcgccttcgccgccccagcgtcaacagtgacctcaagtacatcaaggccaagtccttcgggaaggccaggcccgtagaggggacaagtgcacatccagtgta