

RNF170-ring finger protein 170 Gene

Reference: PTXBC013422

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E3 ubiquitin-protein ligase RNF170; ADSA; SNAX1; ring finger protein 170		

For research use only.

GENE DESCRIPTION	ring finger protein 170
GENE ID	81790
ORIGIN SPECIES	Human
ACCESSIONS	BC013422
VECTOR	pDONR223
NCBI SYMBOL	RNF170
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaaatatacaggtgaagtcaaaagttgaaactggatgatgattcagttatagaaggagtaagcgaccaagtactgtggcagttgtggcagtttgccttgattgctaccctggatatgcacttttcagaaatgtacatcaaaacattcacccagaaaaccaggagctagtaagggtac
ttcgagaacagcttcaaacagacagatgcacctgctgcacctgcagcagttctacactgacatgtactgtcccatctgcctgcaccaagcctcctcccggtggagaccaactgtggacatctttttgtgtaaccttactcctaacagtatttgggtga