

FBXL18-F-box and leucine-rich repeat protein 18 Gene

Reference: PTXBC013435

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Fbl18; F-box/LRR-repeat protein 18; F-box and leucine rich repeat protein 18		
For research use only.			

GENE DESCRIPTION	F-box and leucine-rich repeat protein 18
GENE ID	80028
ORIGIN SPECIES	Human
ACCESSIONS	BC013435
VECTOR	pDONR223
NCBI SYMBOL	FBXL18
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacgcagtcgcccggcgtttggcaagaagtgcgtgtggcggtgcagtcctgtccccccttctcgggcccaggcgtgccccagccctcctcgtgttctgtgtctgtgaagaacctgccctcctggaacacctgcagctgattgggtccaacttctcctccgcatgcccgca
acgagcccgccatccgcaactcgtcccccctgcagccgcgcacagatgtcggggactcggaggtggcccatcggccagctggccttctcgggcacctgacgtcgcacagctgccagcgtccttaagggtccggcgtggtaaatcggcctgcagtcgacagctgtgc
ggcctctcgtcgtggccaactgggcatgatggggaaggtggtgtacatgcccgctctcagacatgttgaagcaactgcaagcggcgtgagggacctcaggctggagcagccctacttcagcggccaacgccagttcttccaggcgtgagccagtgccccctcgtgcagcgcctgtg
cctgggtctcgcagcggcaccctccagcccgatgccgtgctggccttcattggtcgtcgtgcaggtgtgtatgtgtccacctgttcaccggggagtcctcgcacctgcaagagcctgcagcagtcgttctcgcagcttccaggccgagcggcccggttaaacgtcgtcatctcc
ctctcctccagaggcctgaccgacgtcatccgggacgtcccccgtgtgcacctggatga