



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

NTN4-netrin 4 Gene

Reference: PTXBC013591

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRO3091; netrin-4; beta-netrin; hepar-derived netrin-like protein		
For research use only.			

GENE DESCRIPTION	netrin 4
GENE ID	59277
ORIGIN SPECIES	Human
ACCESSIONS	BC013591
VECTOR	pDONR223
NCBI SYMBOL	NTN4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]