

TMEM69-transmembrane protein 69 Gene

Reference: PTXBC013608

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf154; transmembrane protein 69		
For research use only.			

GENE DESCRIPTION	transmembrane protein 69
GENE ID	51249
ORIGIN SPECIES	Human
ACCESSIONS	BC013608
VECTOR	pDONR223
NCBI SYMBOL	TMEM69
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcttgccttcacagaagttttcctcaagcatcttcaagatactgaagtactcttccagtgaggactaagaaccagcagaacagatatacttctctcaagatgtctctccagcaaaactttcccatgtccaaggccttggtttctctcatcttccagcgtatatgagcaagacacagtgct
atcatacatccccctgcagctttaaagcagcagaagcaagcacttctagccagacctcaagcaccatcacttaactgacagcccaagccagcattatgttaactctggcaggactaatcccccttggtgctccaccactggcctgctgacacaaaaacttatatcccatatag
cttttactcagatggccttatggagccagtttctctcttcttgggtggatcagatgggttttctctaccagaaggtagtcagcccaaccagactacctaatttagctagcagtcagctccctcttcttctatggtttgcttctcttattctgaaagacttagtgaagccatagtcacagta
aatgggtatgggagtagcattccacctgaacttttctcttaccacattatcccaactgtttaagccctgaggatagtagtcactttatggccacttttctattataatcacttttagttaaagtagttttccagaaaaaggacataagagacctggtcagataaa