

IL8-interleukin 8 Gene

Reference: PTXBC013615

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IL8; GCP-1; GCP1; LECT; LUCT; LYNAP; MDNCF; MONAP; NAF; NAP-1; NAP1; interleukin-8; T-cell chemotactic factor; alveolar macrophage chemotactic factor I; beta endothelial cell-derived neutrophil activating peptide; beta-thromboglobulin-like protein; chemok		

For research use only.

GENE DESCRIPTION	interleukin 8
GENE ID	3576
ORIGIN SPECIES	Human
ACCESSIONS	BC013615
VECTOR	pDONR223
NCBI SYMBOL	IL8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacttcaaagctggccgtggtctcttggcagccttctgalttctgcagctctgtgtgaagggtcagtttgcgaaggagtgctaaagaacttagatgcagtgcataaagacatactccaaacctttccaccccaaatatcaaagaactgagagtgattgagagtgaccacactgcgcc
aacacagaaattattgtaaagcttctgatggaagagagctctgtctggaccccaaggaaaaactgggtgcagagggttctggagaagtttttgaagagggtgagaattcataa