



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

HXA5-homeobox A5 Gene

Reference: PTXBC013682

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	homeobox protein HOXA5; HOX1; HOX1.3; HOX1C; homeobox protein Hox-A5; homeo box 1C; homeo box A5; homeobox protein Hox-1C; homeobox A5		

For research use only.

GENE DESCRIPTION	homeobox A5
GENE ID	3202
ORIGIN SPECIES	Human
ACCESSIONS	BC013682
VECTOR	pDONR223
NCBI SYMBOL	HOXA5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]