

RNF114-ring finger protein 114 Gene

Reference: PTXBC013695

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E3 ubiquitin-protein ligase RNF114; PSORS12; ZNF313; zinc finger protein 228; zinc finger protein 313; ring finger protein 114		

For research use only.

GENE DESCRIPTION	ring finger protein 114
GENE ID	55905
ORIGIN SPECIES	Human
ACCESSIONS	BC013695
VECTOR	pDONR223
NCBI SYMBOL	RNF114
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgcaacagcgggactgcgggggtgctgcgcagctggcgggggcgggcgaggctgacccctaggacgcttcacgtgtcccggtgtcttagaggtgtacgagaagccggtacaggtgccctgcggacacgtctttgctctgcatgcctgcaggaaatgtctgaagccg
aagaagcctgtctgtgggggtgtgtcgcagcgtctgtgacacctggcgtccgagcgtggagctcgcagcggcgagatcgagagcacagagactcttgcctggtgctgtaagaattctctctgtccaagatccgggtccacgtggctactgttccaaataccagaattacatcatggaagg
tgtgaaggccaccattaaagatgcatctctcagccaaggaatgttccaaacgtttacaccttctgtcttactgtcctgagaagaactttgatcagggaagactgttggaaactgcaaattattccatagcacggataccaaatctgtgtttgtccgatatgtgcctcgatgccctgggga
gaccccaactaccgcagcgccaacttcagagacacatccagcgccggcaccgggtttcttatgacactttgtgattatgatgttgatgaagaggacatgatgaatcaggtgttgacgctccatcatcgaccagtga