



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ADORA2A-adenosine A2a receptor Gene

Reference: PTXBC013780

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

A2aR; ADORA2; RDC8; adenosine receptor A2a; adenosine receptor subtype A2a; adenosine A2a receptor

For research use only.

GENE DESCRIPTION	adenosine A2a receptor
GENE ID	135
ORIGIN SPECIES	Human
ACCESSIONS	BC013780
VECTOR	pDONR223
NCBI SYMBOL	ADORA2A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr