

NHLH1-nescient helix loop helix 1 Gene

Reference: PTXBC013789

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NSCL; NSCL1; bHLHa35; helix-loop-helix protein 1; HEN-1; NSCL-1; class A basic helix-loop-helix protein 35; nescient helix-loop-helix 1		

For research use only.

GENE DESCRIPTION	nescient helix loop helix 1
GENE ID	4807
ORIGIN SPECIES	Human
ACCESSIONS	BC013789
VECTOR	pDONR223
NCBI SYMBOL	NHLH1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgtcaactcagacaccatggagctggacctgcccccacccactcagagactgagtcgggcttcaagtactgtggggcgccggcgccctgaggtgccgggctgggggtccggagggggccaggcccgaggccagagccgggagagcctggccggaagacct
gcagcatctgagccgagggagcgcggcgccggcgccgacagccaagtaccgcacggcccgccacgcgagaaacgcatccgctggaagcctcaacctggccttcgcccagctgcgcaagctgctgctacgtgcccccgacaagaagctctccaagattgagatt
ctgcgctggccatctgtatatctcctacctgaaccagtgctggacgtctga