

PCP4-Purkinje cell protein 4 Gene

Reference: PTXBC013791

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PEP-19; Purkinje cell protein 4; brain specific polypeptide PEP19; brain-specific antigen PCP-4; brain-specific polypeptide PEP-19		

For research use only.

GENE DESCRIPTION	Purkinje cell protein 4
GENE ID	5121
ORIGIN SPECIES	Human
ACCESSIONS	BC013791
VECTOR	pDONR223
NCBI SYMBOL	PCP4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtgagcgacaaaggtgctgggcaaccaatggaaaagacaagacatctggtgaaatgatggacagaagaagtcaagaagaattgacattgacatggatgcaccagagacagaacgtgcagcgggtggccattcagtcagtcagaaaattccagaagaagaaggctgggtctcagtcctag