

BBS10-Bardet-Biedl syndrome 10 Gene

Reference: PTXBC013795

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C12orf58; Bardet-Biedl syndrome 10 protein; Bardet-Biedl syndrome 10		
For research use only.			

GENE DESCRIPTION	Bardet-Biedl syndrome 10
GENE ID	79738
ORIGIN SPECIES	Human
ACCESSIONS	BC013795
VECTOR	pDONR223
NCBI SYMBOL	BBS10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtaccagtgaagttacgaatgggtacttccagagttaccttctcatctatgccagctggtgtgtttgccagtaggtggttaattttgagatctgttacattactatcttcaattatgccaaaaatgccatcaatcagaagaacctgggttagtataatagctaatgcac
ttttaggcattcccaagtcctttataatctaaaacaggaaagtacagctttccacatacatatataagagctgtccatgcactgcaaaccaatcaacccttggaagcagtcagacaggtttggaatcagtaattggtaataaccagctactaactcagttcttcagtttgacaaaaatattaa
ccattgacatggtaatcactgttaagagacaccctcagaagttcacaatcaagattcagaagatgaactataa