

RPSA-ribosomal protein SA Gene

Reference: PTXBC013827

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	37LRP; 67LR; ICAS; LAMBR; LAMR1; LBP; LBP/p40; LRP; LRP/LR; NEM/1CHD4; lamR; p40; 40S ribosomal protein SA; 37 kDa laminin receptor; 37/67 kDa laminin receptor; 67 kDa laminin receptor; colon carcinoma laminin-binding protein; laminin receptor 1 (67kD, ri		

For research use only.

GENE DESCRIPTION	ribosomal protein SA
GENE ID	3921
ORIGIN SPECIES	Human
ACCESSIONS	BC013827
VECTOR	pDONR223
NCBI SYMBOL	RPSA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccggagcccttgatgtcctgcaaatgaaggaggaggtgctccttaagttccttcagcagggaacccacttaggtggccaatctgactccagatggaacagtacatctataaaaggaaaagtgtggcatctatataaaatctcaaggaggacctgggagaagcttctgtgga
gctcgtgcaattgttgcattgaaaacctgctgagtgatgtatctcctcaggaaactggccagaggcgtgctgaagtttgcgtgcccactggagccactccaaatgctggccgcttactctggaacctcactaaccagatccaggcagcctccgggagccacggcttctgt
ggttactgacccaggcgtgaccaccagcctctcacggagcatttatgttaacctacccattgcgctgtgaacacagattcctctgctgtatgtgacattgccatccatgcaacaacaggaggctcactcagtggtttaatgtgttgatgctggctcgggaagttctgcga
tgcgtggcaccattccgtgaacaccatgggaggtcatcctgactctacagagatcctgaagagattgaaaaaagagcagcgtgctgctgagaaggcagtgaccaaggagggaatttcagggtgaatggagctgtcccgtcctgagttcactgctactcagctgaggtt
gcagactggtctgaagggtgacaggtgcccctctgtcctattcagcaattccctactgaagactggagcgtcagcctgccacgggaagactggtctgcagctcccactgctcaggccactgaatgggtaggagcaaccactgactggtcttaa