



Reference: PTXBC013899

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	serine/threonine-protein kinase PLK3; CNK; FNK; PLK-3; PRK; FGF-inducible kinase; cytokine-inducible serine/threonine-protein kinase; proliferation-related kinase; polo like kinase 3		

For research use only.

GENE DESCRIPTION	polo-like kinase 3 (Drosophila)
GENE ID	1263
ORIGIN SPECIES	Human
ACCESSIONS	BC013899
VECTOR	pDONR223
NCBI SYMBOL	PLK3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]