

CCDC124-coiled-coil domain containing 124 Gene

Reference: PTXBC013949

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	coiled-coil domain-containing protein 124; coiled-coil domain containing 124		
For research use only.			

GENE DESCRIPTION	coiled-coil domain containing 124
GENE ID	115098
ORIGIN SPECIES	Human
ACCESSIONS	BC013949
VECTOR	pDONR223
NCBI SYMBOL	CCDC124
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcccaagaagttccaggtgagaacacacagtcggcagcggccggggcagggagggcagagggcagggcggcgtgatgccaagaagcagaaggagctggaggatgcctactggaaggacgacgacaaacacgtcatgaggaggagcagcgaaggaggagaagga
gaagcggcgcctcgaccagctggaacgtaagaaggagacgcagcgcctactggaggaggaggactccaagctcaaggcggcgaaggcggcgggtggccacgtccagcaagggtcaccggggccagatcgaggacacgtgcgccgagaccatcagctcaggaggcccc
ggacacagccgagaaagccaagagccatctggaggctggcgtggaggagaacgtgaaccggcggcgtgctggaggaggcagcgtggaggcgcgacccatcgaggacgccaattgcagtgctcagcgtggcggaggaggcggccgaccggcaccagaagacgcatgagg
cagccttcacagcctttgagggaagccagctgccggctcaacaagagaacccaacatgcggcgtgctgcagctgaaacagctgctcaagaaggagtggtccgctctctgacaacccatgaaccagcggcgctgcccattgcccccaagtga