



Reference: PTXBC013955

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FGFR5; FHFR; fibroblast growth factor receptor-like 1; FGF homologous factor receptor; FGF receptor-like protein 1; FGFR-5; FGFR-like protein; fibroblast growth factor receptor 5		

For research use only.

GENE DESCRIPTION	fibroblast growth factor receptor-like 1
GENE ID	53834
ORIGIN SPECIES	Human
ACCESSIONS	BC013955
VECTOR	pDONR223
NCBI SYMBOL	FGFRL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr