

**UBE2J1-ubiquitin-conjugating enzyme E2, J1 (UBC6 homolog, yeast) Gene**

Reference: PTXBC013973

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGI-76; HSPC153; HSPC205; HSU93243; NCUBE-1; NCUBE1; UBC6; UBC6E; Ubc6p; ubiquitin-conjugating enzyme E2 J1; E2 ubiquitin-conjugating enzyme J1; HSUBC6e; non-canonical ubiquitin-conjugating enzyme 1; ubiquitin-conjugating enzyme E2, J1 (UBC6 homolog, yeas		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	ubiquitin-conjugating enzyme E2, J1 (UBC6 homolog, yeast)
<b>GENE ID</b>	51465
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC013973
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	UBE2J1
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon )**

atggagaccgctacaacctgaaggtccgctgttaaacglttaatgaagaagcggcagaaltgaaagatccaacagatcattaccatgcccagcctttagagataaccttttgaatggcacttcacggttagagggccccagactccgatttgatggaggatttatcacggcgg  
atagtaactccaccagatfaccatgaaccaccaagcattatctcctaaccgctaafggtcgatttgaagtggcagaataatctgtttgagcatcagggccatcctgaaactggcagcctctgtggagtataaggacagcattatagccatcattgggttfaigccaacaaaagga  
gaggagccatagggtctctagattacactctgaggaagaagagcacttgccaaaaatcacaagattctgttgaaggatgtgctctgcatgaaggatgtcctgttgcctttaaacttggaagcgattcaagccaagctgaccaagaagcacaagaactggctagcacaataag  
ctttaaagcagaagtcaattcatctgaaagactatctgagtcagacttaaacactcttttcaactaactgattacaagatgataacctcaaacattccagggtgctacggccagtagatctacggagtcagaaltcctcagcagcatccttcatcaacctaccaactgtagctaaga  
atacctcatgagcctcgacagcgcgggcccagcagcagatcagagaaggtgtctacttcaacagatgtaatccaggccaccagccaagagacaaccaactgatcatggtggcagctgtactgattgtcatcctgacttggcattggcagctcttatattccgacgaatatact  
ggcaaacgaatacatattgactttgagttataa