

WDR92-WD repeat domain 92 Gene

Reference: PTXBC014022

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WD repeat-containing protein 92; WD repeat-containing protein Monad; monad; testicular secretory protein Li 67; WD repeat domain 92		

For research use only.

GENE DESCRIPTION	WD repeat domain 92
GENE ID	116143
ORIGIN SPECIES	Human
ACCESSIONS	BC014022
VECTOR	pDONR223
NCBI SYMBOL	WDR92
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcggccttcgagaagcctcagatcgcgccatccagaaggcctcaactacacgggtgttgactgtaagtgggtgccctgcagcgccaaattgtgacatgggcaacttcgacggggcaccggcgctcattcagctgtacgagatccagcaggggacgtgaagctgttcggg
agattgaaaaggccaaacctataaatgtggaacattgtgcaacatctttacagcagagatatttagctactggagattttgtggaaccttcataatggaatttagaagctccagagatgccagtattctgtaagggccataaagaaattataatgccatagatggcataggtggact
aggaattggagaaggagcacctgaaattgtgactggcagccgagatgggaactgtgaagggtgtgggacccaaggcgaagatgatcctgttgctaataatggaacctgtacaaggagaaaacaagagagactgttggaactgtggcatttggcaatgctataatcaagaagacgtgtgttt
gtgctggctatgacaatggggatcaaaactattgatctcagaatattggcattacgggtgggagacaacatcaaaatgggggtgtgtagcttggagtttgacagaaaagacataagatgaataagtttagtgcacatctctggaaggaaagttccatgttttgacatgagaacacagcat
ccaaccaaaggttttgcctctgtttcagaaaaggctcataatctactgtgtggcaggtccgacacctgccgcagaacaggagagctctttctgacagctggaggcgccggcgccctcacctctggaagtag