



Reference: PTXBC014144

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HUC; HUCL; PLE21; ELAV-like protein 3; ELAV (embryonic lethal, abnormal vision, Drosophila)-like 3 (Hu antigen C); ELAV like neuron-specific RNA binding protein 3; Hu antigen C; hu-antigen C; paraneoplastic cerebellar degeneration-associated antigen; para		

For research use only.

GENE DESCRIPTION	ELAV (embryonic lethal, abnormal vision, Drosophila)-like 3 (Hu antigen C)
GENE ID	1995
ORIGIN SPECIES	Human
ACCESSIONS	BC014144
VECTOR	pDONR223
NCBI SYMBOL	ELAVL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr