



Reference: PTXBC014146

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BETA-4; CPSQ5; SPG47; AP-4 complex subunit beta-1; AP-4 adaptor complex subunit beta; beta 4 subunit of AP-4; beta4-adaptin; spastic paraplegia 47; adaptor related protein complex 4 beta 1 subunit		

For research use only.

GENE DESCRIPTION	adaptor-related protein complex 4, beta 1 subunit
GENE ID	10717
ORIGIN SPECIES	Human
ACCESSIONS	BC014146
VECTOR	pDONR223
NCBI SYMBOL	AP4B1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]