



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

RPL13-ribosomal protein L13 Gene

Reference: PTXBC014167

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

BBC1; D16S444E; D16S44E; L13; 60S ribosomal protein L13; OK/SW-cl.46; breast basic conserved protein 1; ribosomal protein L13

For research use only.

GENE DESCRIPTION	ribosomal protein L13
GENE ID	6137
ORIGIN SPECIES	Human
ACCESSIONS	BC014167
VECTOR	pDONR223
NCBI SYMBOL	RPL13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcggccagccgaatggcgtgcttgaagcccccattccacaagagctgcgcagcggcgctggccacagcggcgccgtaagatccgcagacgtgaagccggcgaagccaagcgcgcgcctatgccccgcgcccgctgggctccatccggcccatc
gtgcgctgccccacggttgggtaccacacagaaggtgctgcgcggcgccgcttcagcttggaggagctcaggggtgcccgcattccacaagaaggtggcccgaccatcggcatttctgtggatccgaggaggcggaacaagctccacggagtcctgcaggccaacgtgcagcggc
cgaagcagatgccgttcccaactcatctcttccccgaagccctgcgcggcccccagaagggagagacttctgctgaagaactgaaactggccaccacgtgacccgcgcgcctatgccgctcgcgcgcctatgccgctcgcgaacgtctatgaagaaggagaagctcgtatctatcactgaggagaagagaagaattt
cgaagcctctgctatcttccgtatgctccctgtcccaacgtcgcgtcttctgcacataagcgcacaaaagcgaaggaagaagccgcgaagaacgtgattgtgaagaagaagaataat