

## PIK3R2-phosphoinositide-3-kinase, regulatory subunit 2 (beta) Gene

Reference: PTXBC014170

|                |                                                                                                                                                                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | MPPH; MPPH1; P85B; p85; p85-BETA; phosphatidylinositol 3-kinase regulatory subunit beta; PI3-kinase subunit p85-beta; PI3K regulatory subunit beta; phosphatidylinositol 3-kinase 85 kDa regulatory subunit beta; phosphatidylinositol 3-kinase, regulatory sub |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | phosphoinositide-3-kinase, regulatory subunit 2 (beta)                                            |
| <b>GENE ID</b>          | 5296                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC014170                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | PIK3R2                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgaagcgtactgcaattgagccttcaatgagactatcaagatcttgaagagcagggccagactcaagagaaatgcagcaaggaatacctggagcgttccggcgtgagggcaacgagaaagagatgcaaggatcctgctgaactccgagcggctcaagtccgcattgccgaga  
tccatgagagccgcacgaagctggagcagcagctgcgggccaggccctggacaacagagagatcgacaagcgcagatgaacagcctcaagccggacctcatgcagctgcgcaagatccgagaccagtacctcgtgtggctcaccagaaaggcccgccgagagaaaaatcaacg  
agtggctgggattaaaaatgagactgaggaccagtacgactcatgaggagcagcagatctcccgaccacgaggaacgcacttggtagctgggcaagatcaaccgcagcagcagagagagatgctgagcggcaagcgggatggcacttctcatccgcgagagcagcca  
gcggggctgctacgctgctcgtggtgtagtgagcggcgacaccaagcactgcgtcatccgcagcgccaccgcttcggcttcgggagccctacaacctgtacgggtcgtgaaggagctggtgctgactaccagcagcctcgtggtgcagcacaacgacgcgtcacgt  
cacctggcgacccagtgcgcgccccgggccccggcccgccgctgcgccccgtga