



Reference: PTXBC014194

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	T-box transcription factor TBX22; ABERS; CLPA; CPX; TBXX; dJ795G23.1; T-box protein 22; T-box 22		

For research use only.

GENE DESCRIPTION	T-box 22
GENE ID	50945
ORIGIN SPECIES	Human
ACCESSIONS	BC014194
VECTOR	pDONR223
NCBI SYMBOL	TBX22
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr