



Reference: PTXBC014195

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C18orf22; HsT169; ribosome binding factor A (putative)		
For research use only.			

GENE DESCRIPTION	chromosome 18 open reading frame 22
GENE ID	79863
ORIGIN SPECIES	Human
ACCESSIONS	BC014195
VECTOR	pDONR223
NCBI SYMBOL	C18orf22
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]