

HOPX-HOP homeobox Gene

Reference: PTXBC014225

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CAMEO; HOD; HOP; LAGY; NECC1; OB1; SMAP31; TOTO; homeodomain-only protein; lung cancer-associated Y protein; not expressed in choriocarcinoma clone 1; odd homeobox protein 1; HOP homeobox		

For research use only.

GENE DESCRIPTION	HOP homeobox
GENE ID	84525
ORIGIN SPECIES	Human
ACCESSIONS	BC014225
VECTOR	pDONR223
NCBI SYMBOL	HOPX
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcggcgagaccgcgagcgcccccacagaggaccaggtggaatcctggagtacaactcaacaaggctgacaaagcccgattccaccacgctgtgcctcatcgcgccgagcgagccttccgaggaggagaccagaaatggttaagcagcgccctggcaaatggc
ggcgtcagaaggcctgcctcagagtcagatccgtcatagactaa