

LOC494150-prohibitin pseudogene Gene

Reference: PTXBC014228

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	prohibitin pseudogene
GENE ID	494150
ORIGIN SPECIES	Human
ACCESSIONS	BC014228
VECTOR	pDONR223
NCBI SYMBOL	LOC494150
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtacagaaccaattatcttgccttctccaccacgtatgtaccaatcatcactggtagcaagatttacagaatgtcaatcacactgcgcacatcttcagcctgttgctagccagcttctcgcacatcttcaccagcatcggagaggactatgatgagcctgtgctgacgtacatcac
gaccgagatcctcaagtcagtggtgctcgtcttgatgctggagaagtattactcagagagagctggtctccaggcaggtgagcaacgaccttacggagcaagcagccacattgggtcatctcggagcagctgtccttgacatactgaccttggaaaggagttcacagaagcagtg
gaagcacaacaggtggtcagcagggaagcagagagggccagatttgaaggaaaaggctgagcagcagaaaaaggctgagcagcagaaaaaggctgagcagcagaaaaaggcagccgtgatctctgctgagggcgactccaaggcaaccgagctgattgtcaactcactggcc
accgcgggggacggcctgatggagctgtgcaagttggaagccgcggaggctctcggaacatga