



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

RHOH-ras homolog gene family, member H Gene

Reference: PTXBC014261

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	rho-related GTP-binding protein RhoH; ARHH; TTF; GTP-binding protein TTF; TTF, translocation three four; ras homolog gene family, member H; ras homolog family member H		

For research use only.

GENE DESCRIPTION	ras homolog gene family, member H
GENE ID	399
ORIGIN SPECIES	Human
ACCESSIONS	BC014261
VECTOR	pDONR223
NCBI SYMBOL	RHOH
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1