

TMEM27-transmembrane protein 27 Gene

Reference: PTXBC014317

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NX-17; NX17; collectrin; kidney-specific membrane protein; transmembrane protein 27		
For research use only.			

GENE DESCRIPTION	transmembrane protein 27
GENE ID	57393
ORIGIN SPECIES	Human
ACCESSIONS	BC014317
VECTOR	pDONR223
NCBI SYMBOL	TMEM27
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttggtgctctctttttctggtgactgccattcatgctgaactctgtcaaccaggtgcagaaaatgcttttaagtgaacttagtatacagaacagctctgggagataaagcatatgcctgggataccaatgaagaataacctctcaagcgatgtagctttctccatgagaaaagttcccaa
cagagaagcaacagaaaattcccatgctcactttgcaatgtaaccagaggggtatcattctggtttgtggttacagacccttcaaaaatcacaccctctcgtctgtgaggtgcaatcagccataagaatgaacaagaaccggatcaacaatgcctctttctaaatgaccaaactctggaatttt
taaaaatcccttccacttgaccacccatggaccatctgtgccatctggattattatatttgggtgatattttgcatcatcatagttgcaattgcactactgattttatcagggaatctggcaacgtagaagaagaacaagaaccatctgaagtgatgacgctgaagataagttgaaaac
atgatcacaattgaaaatggcatccctctgatccctggacatgaaggagggcatattaatgatgccttcacacagaggatgagagggtcacccctctctga