

STMN1-stathmin 1/oncoprotein 18 Gene

Reference: PTXBC014353

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf215; LAP18; Lag; OP18; PP17; PP19; PR22; SMN; stathmin; leukemia-associated phosphoprotein p18; metablastin; oncoprotein 18; phosphoprotein 19; phosphoprotein p19; prosolin; stathmin 1/oncoprotein 18; testicular tissue protein Li 189; transmembrane p		

For research use only.

GENE DESCRIPTION	stathmin 1/oncoprotein 18
GENE ID	3925
ORIGIN SPECIES	Human
ACCESSIONS	BC014353
VECTOR	pDONR223
NCBI SYMBOL	STMN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcttctgatatccaggtgaaagaactggagaagcgtgcctcaggcaggttttgagctgattctcagccctcggcctcagaaagaatctgtccagaattcccccttccccctcagaagaagagattttccctggaggaaattcagaagaattagaagctgcagaagaagacgca
agtcctcatgaagctgaggtcttgaagcagctggctgagaacgagagcagcagagaagaagtgcttcagaaggcaatagaagaacaacaacttcagtaaatggcagaagagaactgacccacaaatggaagctaatagaagaaccgagaggcacaatggctgccaactg
gaacgtttgcgagagaaggataagcacattgaagaagtcggaagaacaagaatccaagacccctgctgacgagactgaagctgactaa