

OR7E91P-olfactory receptor, family 7, subfamily E, member 91 pseudogene Gene

Reference: PTXBC014374

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	olfactory receptor, family 7, subfamily E, member 91 pseudogene
GENE ID	79315
ORIGIN SPECIES	Human
ACCESSIONS	BC014374
VECTOR	pDONR223
NCBI SYMBOL	OR7E91P
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtcttctctccaaacctgtgtctgggtgacatcggtttcacctggccacgggtcccaagatgattgtggacatcggttcgcatagcggagtcacatcttatcggaactgcctgacacggatgtcttcttggctcttttgcattgtgtagaagacatgctcctgactgtgatggcctatgac
tggtttgtagccatctgtcgccctctgacacccagtcacgtgaatctcacctctgtgtctcttagtttcgggtgctcttttcccttagcctgttggtatccacgtgcgcagttggattgtgttgcaattcacctcttcaagaatgtggaatctctaatttctgtgacccatctcaacctctcaagct
tgccctgtctgacagcatcatgtagcatgttcataatttcgtagtactatgtttgttttcccaatttcaggatcctttgtcttactataaaattgtccctccattctaa