

SURF2-surfeit 2 Gene

Reference: PTXBC014411

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SURF-2; surfeit locus protein 2; surfeit 2		
For research use only.			

GENE DESCRIPTION	surfeit 2
GENE ID	6835
ORIGIN SPECIES	Human
ACCESSIONS	BC014411
VECTOR	pDONR223
NCBI SYMBOL	SURF2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcaggttgcggcgacgtgcggcggtttctgcgggagcacccgagcctgcggctccagacggagcccgcaagggtgaggtgcacctgcacaggtcacgagctgccctgccgctgccgagctccaggtctacacccgcgcaaaaagttaccagcggctggtccgcgct
ccccggccttcgactatgcagagttcgcggcgacatcgtgccagcaccagaacccgacagttgttctgcaaaactacccctgcggcaccatcaacaagtgcacagaacacgtgctgaggcacacccaggccggcggtaccagcagctctgttaaatatgaagaatgtcagaa
gcaaggggtggagtacgtgcctgcctgcctggtgcacgggagggagggagggagggaccagatggacgtgacggcctgcggcggaagccttctggagcccacatccagtgtgagggggagctgcaagtgatgacagcatgacagacctgtaccacctgagctatc
accagaaaggaccttggaagcacggaggtgggatggcactgatgactttttgacagacaaaggatgagaaggcaagcccccaagagagaaggccactgatgaggcgaggagagacgacctgtaccgagggctggtccagaagcggggaagaagcagttgggctc
gttgaataaagaagttcaagagtcacccgcaaaccaagagcttcagctcctgtaaacagccaggttaa