



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

STK10-serine/threonine kinase 10 Gene

Reference: PTXBC014413

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRO2729; serine/threonine-protein kinase 10; lymphocyte-oriented kinase; serine/threonine kinase 10		

For research use only.

GENE DESCRIPTION	serine/threonine kinase 10
GENE ID	6793
ORIGIN SPECIES	Human
ACCESSIONS	BC014413
VECTOR	pDONR223
NCBI SYMBOL	STK10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcattctctctgcttctccagctttgctttccacttaagacaagacgttctctatgtgggggggggggacgggggaaagaggagctgaattttattctctcttcccgcttttccatctcgcgtccccctctctgcacatgggtgtaatcacacatacgcgtacacacaggaactctctcgtggcgacctctctgcacagctgggttggaacacgttctctgcagctctgcccactgagcatctttatgttgagacagctgtgtgtgcacgtctcgggtttggcaccggagcgctgaaagacagagctgtaactgccagctcgatcgtctgctgcccccacagcgtctctctgacagcagctgattctctcacacacagcgtctccacacccactctctccatag