

RPL23A-ribosomal protein L23a Gene

Reference: PTXBC014459

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L23A; MDA20; 60S ribosomal protein L23a; melanoma differentiation-associated gene 20; ribosomal protein L23a		

For research use only.

GENE DESCRIPTION	ribosomal protein L23a
GENE ID	6147
ORIGIN SPECIES	Human
ACCESSIONS	BC014459
VECTOR	pDONR223
NCBI SYMBOL	RPL23A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccgaaagcgaagaaggaagctcctgccctcctaaagctgaagccaaagcgaaggcttaagccaaaggaagcagtggtgaaaggtgtccacagccacaaaaagaagaagatccgcacgtcaccacattccggcgccgaagacactgcgactccggagacagccca
aatatcctcgggaagagcgtcccgagagaacaagctgaccactatgctatcatcaagtttccgctgaccactgagctcgcctgaagaagatagaagacaacaacacactgtgttcattgtggtgttaagccaacaagcaccagattaaacaggctgtgaagaagctgtatgacattg
atgtggccaaggtaacaccctgattcggcctgatggagagaagaaggcatatgttcgactggctcctgattacgagcttggatgttgccaacaaattgggatcatctaa