

SEC11A-SEC11 homolog A (*S. cerevisiae*) Gene

Reference: PTXBC014508

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	signal peptidase complex catalytic subunit SEC11A; 1810012E07Rik; SEC11L1; SPC18; SPCS4A; sid2895; SEC11 homolog A; SEC11-like protein 1; SPase 18 kDa subunit; endopeptidase SP18; microsomal signal peptidase 18 kDa subunit; signal peptidase complex (18kD)		

For research use only.

GENE DESCRIPTION	SEC11 homolog A (<i>S. cerevisiae</i>)
GENE ID	23478
ORIGIN SPECIES	Human
ACCESSIONS	BC014508
VECTOR	pDONR223
NCBI SYMBOL	SEC11A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgtcttagacttttggacgatgtgcggcgatgaacaagcgagccttattatcaagtcctaaatttggatgattgtctcatcgccactaatgatctggaaggggttaaggtaataactggaagtgaagtcgattgtagtgtgctcagtgccagcatggaacctgcattcata
gaggagatcttcttctaacaatcgagttgaagatccatagtgaggagaaattgtgttttaggatagaaggaagagagattcctatagttcaccgagtcctgaagattcatgaaagcaaaatgggcataatcaagttttgaccaaaggagataataatgcggttgatgaccgaggcc
tctataacaaggacaacattggctagagaaaaagatgttggggagagccagggaattgttccttatattggaattgtgacgatcctcatgaatgactatcctaaatttaagtatgcagttctcttttctgggtttattcgtgctggttcacgtgagtaa