



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.ProteoGenix.science

ACRV1-acrosomal vesicle protein 1 Gene

Reference: PTXBC014588

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	D11S4365; SP-10; SPACA2; acrosomal protein SP-10; sperm protein 10; acrosomal vesicle protein 1		

For research use only.

GENE DESCRIPTION	acrosomal vesicle protein 1
GENE ID	56
ORIGIN SPECIES	Human
ACCESSIONS	BC014588
VECTOR	pDONR223
NCBI SYMBOL	ACRV1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1