

C20orf141-chromosome 20 open reading frame 141 Gene

Reference: PTXBC014591

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C20orf141; dJ860F19.4; chromosome 20 open reading frame 141		

For research use only.

GENE DESCRIPTION	chromosome 20 open reading frame 141
GENE ID	128653
ORIGIN SPECIES	Human
ACCESSIONS	BC014591
VECTOR	pDONR223
NCBI SYMBOL	C20orf141
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacccggctctgttaccagaccgaagcagctgaggatccgatccagttctccaaggggctgggtgctggggaggggtcaggtagtcagtcgctccactgtatccactggggccctagctgggccagctcctggacagtgtcctatggctggggcactaggactgac
aatccaggcagctctttccaccactggccagccctgctgctgtctgtcagcttctcctcactttgacctgctccatagggccgaggtcacactctccacagcgcgaactctccaccagggccagagtcagggggccggtgaaggtcctggacagcaggaggctactcctgcaa
atgggtacagctcaggacaacttagcctccaggacgactgctgctgctcatggggctggggccgctcctgagagcctgtggcatgccctgacccctgttggcctggctttctgctccatccttgggcctga