

FAM92A1-family with sequence similarity 92, member A1 Gene

Reference: PTXBC014598

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM92A1; FAM92A1; family with sequence similarity 92 member A1; family with sequence similarity 92 member A		

For research use only.

GENE DESCRIPTION	family with sequence similarity 92, member A1
GENE ID	137392
ORIGIN SPECIES	Human
ACCESSIONS	BC014598
VECTOR	pDONR223
NCBI SYMBOL	FAM92A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgagcgcaccctgaaaaccggaacgctcaaacgaaacaactgcaaacagctgtcctcaaatgtggagaagcatttggagaactgtgccaatctcgtgctctatgtcggaaaactgccaggctgagagacaaagcagacctctggtgaatgaaattaacgcgtatctgcta
cagagaccccgcatftaaagctggcctgatgaacttgcagatgagttgccaactcagattatcacaagcagagggtgaaagactgaaagcaaaagtagtgaaccctgaaaactatggaccattgtgaaaatgaaacgggatgacctaaagcaacactcacagcaaggaaat
cgagaagctaagcaattaactcagttagaagaacacgtcagcgaacccatctgatcgacatgttatttcacagggcagaacggaattacagagagctgcaatggatgctagccgaacaagctgcatctggagaaactatcaacttgaaggcagaaaatgaaggatataaaga
ctatatttctgaatttatacaatcgaatgttatttcacggcaagcgtttagaggtctacactgctgcctaccagaatacaaaaacattgatgaagatgaagatttagaggtttccgaaattctgtatgcaccagattatcctcgtttagatattgtaagacaaatcaaagtcaccttca
gagatcactgctcagctaaagtgtatctggaacaggacaggaggctgctggccgcaactgccaacccctgatctag