

C1orf182-chromosome 1 open reading frame 182 Gene

Reference: PTXBC014605

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf182; SIP; SSTK-IP; TSSK6-activating co-chaperone protein; SSTK-interacting protein (SSTK-IP); TSSK6 activating co-chaperone; TSSK6 activating cochaperone		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 182
GENE ID	128229
ORIGIN SPECIES	Human
ACCESSIONS	BC014605
VECTOR	pDONR223
NCBI SYMBOL	C1orf182
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagcggcacactagtcctaacagaaaaagttccagccaaagagggaagctaagctgtgcctctctgtagagcaaaacctccccagctatattaatctcaagcaagttccccaccagccactttctgaacatccagacaacaaagctgcctcggttgatcacagcccaaggaa
tgctaggactcctggaatgtatgtatgcaaacctccagctcagaccagctcgcaccaacagatggctgtttggaacatttacaggcatctgtgacacaactggctcctggaggagggaagcaataactctctctccagccttatctcctaatacattgttaatacctgcccgaattca
gtaaatga