

HMGN2-high-mobility group nucleosomal binding domain 2 Gene

Reference: PTXBC014644

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HMG17; non-histone chromosomal protein HMG-17; high mobility group nucleosome-binding domain-containing protein 2; high mobility group protein N2; high-mobility group (nonhistone chromosomal) protein 17; nonhistone chromosomal protein HMG-17; high mobility		

For research use only.

GENE DESCRIPTION	high-mobility group nucleosomal binding domain 2
GENE ID	3151
ORIGIN SPECIES	Human
ACCESSIONS	BC014644
VECTOR	pDONR223
NCBI SYMBOL	HMGN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcccaagagaaaggctgaaggggatgtaaggagataaagcaagggtgaaggacgaaccacagagaagatccgcgaggttgctgctaaacctgtcctccaaagccagagcccaagcctaaaaagggccctgcaaagaaggagagagaaggtacccaaagggaagggg
aaaagctgatgctggcaggagggaataacccctgcagaaaatggagatgcaaaacagaccaggcacagaaaagctgaaggtgctggagatgcccaagtga