

HOXB6-homeobox B6 Gene

Reference: PTXBC014651

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HOX2; HOX2B; HU-2; Hox-2.2; homeobox protein Hox-B6; homeo box 2B; homeo box B6; homeobox protein Hox-2.2; homeobox protein Hox-2B; homeobox protein Hu-2; homeobox B6		

For research use only.

GENE DESCRIPTION	homeobox B6
GENE ID	3216
ORIGIN SPECIES	Human
ACCESSIONS	BC014651
VECTOR	pDONR223
NCBI SYMBOL	HOXB6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagttcttattctgaactccacctccccgtcactctgcccagcgggcaggagtcctctctggccagctaccgctctattctgctgggctatgcggaccgctgagacattaccccgccctacgggcccaggcggccaggacaaggcttgcacttctctattaccgcc
ggcgggctgtgctacggccgagcgcgcgcctgcactacggccgcccgcgccttctaccgcgagaaagagtcggcctgcgcactctccggcggcagcagcagcccccgttccaccccgagccgcggaagtcggactgcgcaggaagagcgtgttcggcgagaca
gaagagcagaagtgctccactccggtctacccgtgcatgcagcggatgaattcgtgcaacagttctcttgggcccagcggccgagggccgacacacacgttaccagacgtggagctggagaaggagtttcaactacaatcgtactctgacgcggcggcgcatcg
agatcgcgcacgccctgtcctgacggagaggcagatcaagatatggtccagaaccgacgcatgaagtggaaaaaggagagcaaacgtctcagcgcgtctcagctcagtcggagaggagggaagaaaaacaggccgagtgga