



Reference: PTXBC014654

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ABP-280; ABPX; CSBS; CVD1; FLN; FLN-A; FLN1; FMD; MNS; NHBP; OPD; OPD1; OPD2; XLVD; XMVD; filamin-A; actin binding protein 280; alpha-filamin; endothelial actin-binding protein; filamin A, alpha; filamin-1; non-muscle filamin; filamin A		

For research use only.

GENE DESCRIPTION	filamin A, alpha (actin binding protein 280)
GENE ID	2316
ORIGIN SPECIES	Human
ACCESSIONS	BC014654
VECTOR	pDONR223
NCBI SYMBOL	FLNA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]