

C12orf43-chromosome 12 open reading frame 43 Gene

Reference: PTXBC014661

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C12orf43; chromosome 12 open reading frame 43		
For research use only.			

GENE DESCRIPTION	chromosome 12 open reading frame 43
GENE ID	64897
ORIGIN SPECIES	Human
ACCESSIONS	BC014661
VECTOR	pDONR223
NCBI SYMBOL	C12orf43
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgcccgagtgacgtgagcattcggaagtagtaacagcagtagcgaatcgaggagctggagcgtgcccggaggcggaatgccggcttggggctggagcaacgcccgcacgtggcaggaagccaagagccggtgctgcaaatagccagttgtcaacctcc
caaccgagcctcaggcataagtgatgagcatgaacaagatggcaacgagcttcagaccacccctgaattccgagccacgtagccaagaagctgggagccctgctggacagcttcattaccatctcagaagcagcaaaaggagccagcaaaagtaaggtacagaagtcgcttgg
aggatgatggttccgctttctcacatctgtccctggaggccgtgagaagggaagagctcccccaacccgcccgaagcgacagccctccagctccagcagtgaggacagtgacgaggagtgccggcgtgcccgggagggcagctgtgtcggcgtccgacatcctacaggagtcagc
cattcacgcccctggaaactgtggagaaggaggcaagaagaaagggaagtfgaaaaaagaaagccaagaagtgccagctgtcgtgctgctgccaccacccaccagcatggccacagtcagaagcagaagtcagggtgagctcaacggggaccaggtgtcgttggg
accaaaaagaagaaaggcaagaaggccagcagacctctccattccaccagcaaaagtgctacagctatacctgcaaaactga