



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ZNF524-zinc finger protein 524 Gene

Reference: PTXBC014666

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

zinc finger protein 524

For research use only.

GENE DESCRIPTION	zinc finger protein 524
GENE ID	147807
ORIGIN SPECIES	Human
ACCESSIONS	BC014666
VECTOR	pDONR223
NCBI SYMBOL	ZNF524
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]