



Reference: PTXBC014789

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GALNTL2; GALNACT15; GALNT13; GALNT7; PIH5; pp-GalNAc-T15; polypeptide N-acetylgalactosaminyltransferase 15; UDP-GalNAc transferase T15; UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase-like protein 2; UDP-N-acetyl-alpha-D-galactosamine:polypeptide		

For research use only.

GENE DESCRIPTION	UDP-N-acetyl-alpha-D-galactosamine:polypeptide acetylgalactosaminyltransferase-like 2	N-
GENE ID	117248	
ORIGIN SPECIES	Human	
ACCESSIONS	BC014789	
VECTOR	pDONR223	
NCBI SYMBOL	GALNTL2	
LENGTH	ProteoGenix bp	
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG	

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr