



Reference: PTXBC014844

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	segment polarity protein dishevelled homolog DVL-2; dishevelled 2 (homologous to <i>Drosophila</i> dsh); dishevelled, dsh homolog 2; dishevelled segment polarity protein 2		

For research use only.

GENE DESCRIPTION	dishevelled, dsh homolog 2 (Drosophila)
GENE ID	1856
ORIGIN SPECIES	Human
ACCESSIONS	BC014844
VECTOR	pDONR223
NCBI SYMBOL	DVL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]