

RABL2B-RAB, member of RAS oncogene family-like 2B Gene

Reference: PTXBC014879

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	rab-like protein 2B; RAB, member of RAS oncogene family-like 2B		
For research use only.			

GENE DESCRIPTION	RAB, member of RAS oncogene family-like 2B
GENE ID	11158
ORIGIN SPECIES	Human
ACCESSIONS	BC014879
VECTOR	pDONR223
NCBI SYMBOL	RABL2B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagaagacaaacaaaccgagtgagttggaccaagggaagtatgatgctgatgacaacgtgaagatcatctgcctgggagacagcgcagtgaggcaaatccaaactcatggagagatttctcatggatggcttcagccacagcagctgtccacgtacgcctgaccctgtacaa
gcacacagccacggttagatggaaggaccatccttgggactttgggacacggcagccaggagcgttccagagcatgcatgctctactaccacaaggccacgcctgcatcatggtgtttgatgtacagaggaaagtcacctataggaaacctgagcacctggtatacagagcttcg
ggagttcaggccagagatcccatgcatcgttggccaataaattgatgcagacataaacgtgacccaaaaagcttcaatttggcaagaagtctccctgcccctgtatttcgtctcggctgctgatgtaccaatgttgaagctcttcaatgatgcaattcgattagctgtgtcttacaac
agaactcccaggacttcaggatgagattttcaggagctcgagaactcagcttggagcagggaagaggaggacgtgccagaccagggaacagagcagcagcagcagacccatcagaggaggcgccctctccccacagctga