



Reference: PTXBC014927

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PHLPP; PLEKHE1; PPM3A; PH domain leucine-rich repeat-containing protein phosphatase 1; PH domain-containing family E member 1; SCN circadian oscillatory protein; pleckstrin homology domain containing, family E (with leucine rich repeats) member 1; protein		

For research use only.

GENE DESCRIPTION	PH domain and leucine rich repeat protein phosphatase
GENE ID	23239
ORIGIN SPECIES	Human
ACCESSIONS	BC014927
VECTOR	pDONR223
NCBI SYMBOL	PHLPP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]